


Product Datasheet

TNF alpha Recombinant Protein (HEK293 Expressed) (orb1227854)

Description	TNF alpha Recombinant Protein (HEK293 Expressed)
Tested Applications	WB
Target	TNF
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Val 77 - Leu 233
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	17.4 kDa
Uniprot ID	P01375
NCBI	P01375
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human
TNF-alpha,
Tag Free,
low
endotoxin...