


Product Datasheet

SIRPA Recombinant Protein (orb1227872)

Description	SIRPA Recombinant Protein
Tested Applications	WB
Target	SIRPA
Preservatives	50 mM tris, 100 mM glycine, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	C-mFc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Glu 31 - Arg 370
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	63.5 kDa
Uniprot ID	P78324
NCBI	NP_001035111
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of
rh
CD172a/SIRPA
mouse Fc C...