


Product Datasheet

SIRP alpha Recombinant Protein (orb1227874)

Description	SIRP alpha Recombinant Protein
Tested Applications	WB
Target	SIRPA
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Glu 31 - Arg 370
Purity	> 95% as determined by SDS-PAGE and SEC-HPLC.
Source	HEK293 cells
MW	38.1 kDa
Uniprot ID	P78324
NCBI	NP_001035111
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human SIRP
alpha, His
Tag (HPLC-
verified...