


Product Datasheet

SIGLEC6 Recombinant Protein (orb1227875)

Description	SIGLEC6 Recombinant Protein
Tested Applications	WB
Target	SIGLEC6
Preservatives	50 mM tris, 100 mM glycine, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	C-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Gln 16 - Val 320
Purity	> 92% as determined by SDS-PAGE.
Source	HEK293 cells
MW	60.3 kDa
Uniprot ID	O43699
NCBI	AAH35359
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of rh
SIGLEC6/CD327
Fc Chimera...