


Product Datasheet

Semaphorin 4A Recombinant Protein (orb1227878)

Description	Semaphorin 4A Recombinant Protein
Tested Applications	WB
Target	SEMA4A
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Gly 33 - His 683
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	72.7 kDa
Uniprot ID	Q9H3S1
NCBI	AAH20974
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human
Semaphorin
4A, His Tag
on SDS-
PAGE...