


Product Datasheet

CD124 Recombinant Protein (orb1227978)

Description	CD124 Recombinant Protein
Tested Applications	WB
Target	IL4R
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Met 26 - His 232
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	24.5 kDa
Uniprot ID	P24394
NCBI	NP_000409.1
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



Human IL-4
R alpha, His
Tag on SDS-
PAGE ...