


Product Datasheet

IL18 BP Recombinant Protein (orb1227981)

Description	IL18 BP Recombinant Protein
Tested Applications	WB
Target	IL18BP
Preservatives	Tris with Glycine, Arginine and NaCl, pH 7.5
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	C-Fc Tag
Note	For research use only
Application notes	This recombinant protein can be used for WB. For research use only.
Protein Sequence	Thr 31 - Gly 194
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	44.3 kDa
Uniprot ID	O95998
NCBI	NP_001138529
Dilution Range	This recombinant protein can be used for WB. For research use only.
Expiration Date	6 months from date of receipt.



The purity of
rh IL18BP/IL-
18BPα Fc
Chim...