


Product Datasheet

Uromodulin antibody (orb12281)

Description	Goat polyclonal antibody to Uromodulin
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Uromodulin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 65kDa band observed in lysates of cell line HeLa (calculated MW of 69.8kDa according to NP_003352.2). Recommended concentration: 0.03-0.1µg/ml.
Protein Sequence	PITRKGVQATVS
Clonality	Polyclonal
MW	69.8
Entrez	7369
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using Uromodulin antibody (orb12281).