


Product Datasheet

Uromodulin antibody (orb12282)

Description	Goat polyclonal antibody to Uromodulin
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	Uromodulin (aa525-536)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the Plasma membrane and Golgi of HeLa cells. Recommended concentration: 10µg/ml.
Protein Sequence	DRCPHTRDSTIQ
Clonality	Polyclonal
MW	69.8
Entrez	7369
Dilution Range	ELISA: 1:64000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...