

Product Datasheet

SARS-CoV-2 (COVID-19) Helicase Recombinant Protein (orb1228226)

Description	SARS-CoV-2 (COVID-19) Helicase Recombinant Protein
Target	ORF1ab
Preservatives	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.5.
Form/Appearance	Liquid
Storage	Store below -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	6His tag at the C-terminus
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Ala1-Gln601
Purity	Greater than 80% as determined by reducing SDS-PAGE.
Source	E. coli
Uniprot ID	P0DTC1
NCBI	YP_009725297.1
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.