

Product Datasheet

SARS-CoV-2 (COVID-19) S Protein RBD (N501Y) Recombinant Protein (orb1228229)

Description	SARS-CoV-2 (COVID-19) S Protein RBD (N501Y) Recombinant Protein
Target	S
Preservatives	Supplied as a 0.2 uM filtered solution of PBS, pH 7.4
Form/Appearance	Liquid
Storage	Reconstituted protein solution should be store below -20°C.
Tag	6His tag at the C-terminus
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Arg319-Ser541(N501Y)
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	HEK293 cells
MW	26kDa; 33-36 kDa, reducing conditions
Uniprot ID	QHD43416.1
NCBI	QHD43416.1
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.