


## Product Datasheet

# SARS-CoV-2 (COVID-19) S1 (D614G) Recombinant Protein (orb1228234)

|                          |                                                                                    |
|--------------------------|------------------------------------------------------------------------------------|
| <b>Description</b>       | SARS-CoV-2 (COVID-19) S1 (D614G)<br>Recombinant Protein                            |
| <b>Target</b>            | S                                                                                  |
| <b>Preservatives</b>     | Supplied as a 0.2 uM filtered solution of PBS, pH 7.4.                             |
| <b>Form/Appearance</b>   | Liquid                                                                             |
| <b>Storage</b>           | Reconstituted protein solution should be store below -20°C.                        |
| <b>Tag</b>               | C-10 His tag                                                                       |
| <b>Note</b>              | For research use only                                                              |
| <b>Application notes</b> | This recombinant protein can be used for biological assays. For research use only. |
| <b>Protein Sequence</b>  | Gln14-Arg685(D614G)                                                                |
| <b>Purity</b>            | Greater than 90% as determined by reducing SDS-PAGE.                               |
| <b>Source</b>            | Mammalian Cells                                                                    |
| <b>MW</b>                | 76.6kDa, 100-130kDa                                                                |
| <b>Uniprot ID</b>        | QHD43416.1                                                                         |
| <b>NCBI</b>              | <b>QHD43416.1</b>                                                                  |
| <b>Dilution Range</b>    | This recombinant protein can be used for biological assays. For research use only. |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                     |



Immobilized  
Human  
ACE-2-Fc at  
5 ug/ml  
(1...