


Product Datasheet

SARS-CoV-2 (COVID-19) S Protein RBD-SD1 Recombinant Protein (orb1228237)

Description	SARS-CoV-2 (COVID-19) S Protein RBD-SD1 Recombinant Protein
Target	S
Preservatives	Supplied as a 0.2 uM filtered solution of PBS, pH 7.4.
Form/Appearance	Liquid
Storage	Reconstituted protein solution should be store below -20°C.
Tag	C-Avi-6His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Arg319-Ser591
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Mammalian Cells
MW	33.2kDa, 38-35kDa
Uniprot ID	QHD43416.1
NCBI	QHD43416.1
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.



Immobilized
Human
ACE-2-Fc at
2 ug/ml
(1...