

Product Datasheet

SARS-CoV-2 (COVID-19) S1 Recombinant Protein (orb1228243)

Description	SARS-CoV-2 (COVID-19) S1 Recombinant Protein
Target	S
Preservatives	Supplied as a 0.2 uM filtered solution of PBS, pH 7.4
Form/Appearance	Liquid
Storage	Reconstituted protein solution should be store below -20°C.
Tag	C-mFc tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Gln14-Arg685
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Mammalian Cells
MW	101.3kDa, 130-160kDa
Uniprot ID	QHD43416.1
NCBI	QHD43416.1
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.