


Product Datasheet

SLC26A6 antibody (orb12284)

Description	Goat polyclonal to SLC26A6.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SLC26A6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize isoform 1, 2, 3 and 4 NP_075062.2;NP_599025.2;NP_602298.2;NP_001035544.1. Please note that this antibody was shown unfit for Mouse tissues both in IHC and WB.
Protein Sequence	SAPRTHQWRT
Clonality	Polyclonal
MW	83.0; 82.8; 80.8; 80.9
Entrez	65010
Dilution Range	ELISA: 1:8000, WB: 10 µg/ml, IHC-P: 0.1-0.3mg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human renal cor...