

Product Datasheet

Aminoacylase-3 Recombinant Protein (orb1228643)

Description

Aminoacylase-3 Recombinant Protein

Target	ACY3
Preservatives	Supplied as a 0.2 um filtered solution of 20mM Tris,100mM NaCl,1mM DTT,10%Glycerol,pH 8.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	N-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Met1-Ser319
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	E. coli
MW	37.4 kD
Uniprot ID	Q96HD9
NCBI	Q96HD9
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.