



Product Datasheet

Calbindin antibody (orb12289)


Description	Goat Polyclonal antibody to Calbindin
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	calbindin D28
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx. 26kDa band observed in Human and Rat Kidney lysates and in Rat Brain lysates (calculated MW of 30.0kDa according to Human NP_004920.1 and Rat NP_114190.1). The observed molecular weight is routinely observed by other sources. Recommended concentrat. IHC: Paraffin embedded Human Kidney. Recommended concentration: 5-7µg/ml.
Protein Sequence	KTFVDQYGRDDGK
Clonality	Polyclonal
MW	30
Entrez	793
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Hum...



Immunohistochemical analysis of staining...



Western blot analysis of staining of Mou...