

Product Datasheet

Endoplasmic Reticulum Resident Protein 44 Recombinant Protein (orb1228999)

Description

Endoplasmic Reticulum Resident

Target	ERP44
Preservatives	Supplied as a 0.2 um filtered solution of 20mM Tris HCl 10% glycerol, PH7.5. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	C-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Glu30-Asp402
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Human Cells
MW	44.2 kD
Uniprot ID	Q9BS26
NCBI	Q9BS26
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.