

Product Datasheet

TWEAK Receptor Recombinant Protein (orb1229021)

Description

TWEAK Receptor Recombinant Protein

Target	TNFRSF12A
Preservatives	Lyophilized from a 0.2 um filtered solution of 20mM PB,150mM NaCl,pH7.4. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Glu28-Trp79
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Human Cells
MW	35.6 kD
Uniprot ID	Q9NP84
NCBI	Q9NP84
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.