

Product Datasheet

Sorbitol Dehydrogenase Recombinant Protein (orb1229232)

Description
Sorbitol Dehydrogenase

Target	SORD
Preservatives	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 0.2M NaCl, 5mM DTT, 20% glycerol, pH8.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	C-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Ala2-Pro357
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Human Cells
MW	39.3 kD
Uniprot ID	Q00796
NCBI	Q00796
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.