

Product Datasheet

15-PGDH Recombinant Protein (orb1229241)

Description	15-PGDH Recombinant Protein
Target	HPGD
Preservatives	Supplied as a 0.2 um filtered solution of 20mM HEPES,150mM NaCl,pH7.4. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	C-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Met1-Gln266
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Human Cells
MW	30 kD
Uniprot ID	P15428
NCBI	P15428
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.