



Product Datasheet

ENC1 antibody (orb12293)

Description	Goat polyclonal antibody to ENC1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ENC1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in Rat Ovary lysates (calculated MW of 66.1kDa according to Rat NP_001003401.1. Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Mouse Brain for this particular batch.. IHC: Paraffin embedded Human Adrenal Gland. Recommended concentration: 5µg/ml.
Protein Sequence	TSVSHDKLPKVQC
Clonality	Polyclonal
MW	66.1
Entrez	8507
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb12293
(1µg/ml)
staining of
Rat Ovary...



orb12293
(5µg/ml)
staining of
paraffin ...