


Product Datasheet

MCTS1 antibody (orb12294)

Description	Goat polyclonal to MCTS1.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MCTS1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms NP_054779.1;NP_001131026.1.
Protein Sequence	DEKENVSNCIQLKTS
Clonality	Polyclonal
MW	20.6; 20.6
Entrez	28985
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of MOLT4 lysate us...