

Product Datasheet

CD355 Recombinant Protein (orb1229573)

Description

CD355 Recombinant Protein

Target

CRTAM

Preservatives

Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.

Form/Appearance

Lyophilized

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Tag

C-6 His tag

Note

For research use only

Application notes

This recombinant protein can be used for biological assays. For research use only.

Protein Sequence

Ser18-Ser286

Purity

Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.

Source

Human Cells

MW

30.99 kD

Uniprot ID
O95727
NCBI
O95727
Dilution Range

This recombinant protein can be used for biological assays. For research use only.

Expiration Date

6 months from date of receipt.