

Product Datasheet

Complement Factor B Recombinant Protein (orb1229575)

Description

Complement Factor B Recombinant

Target	CFB
Preservatives	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0. It is not recommended to reconstitute to a concentration less than 100 ug/ml.
Form/Appearance	Liquid
Storage	Store at -20°C, stable for 6 months after receipt. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Tag	C-6 His tag
Note	For research use only
Application notes	This recombinant protein can be used for biological assays. For research use only.
Protein Sequence	Thr26-Leu764
Purity	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
Source	Human Cells
MW	84.07 kD
Uniprot ID	P00751
NCBI	P00751
Dilution Range	This recombinant protein can be used for biological assays. For research use only.
Expiration Date	6 months from date of receipt.