

Product Datasheet

C1qTNF1 Recombinant Protein (orb1229596)

Description

C1qTNF1 Recombinant Protein

Target

C1QTNF1

Preservatives

Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.

Form/Appearance

Lyophilized

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at 4-7°C for 2-7 days.
Aliquots of reconstituted samples are stable at -20°C for 3 months.

Tag

C-6 His tag

Note

For research use only

Application notes

This recombinant protein can be used for biological assays. For research use only.

Protein Sequence

Arg26-Pro281

Purity

Greater than 95% as determined by reducing SDS-PAGE.
Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.

Source

Human Cells

MW

30.24 kD

Uniprot ID
[Q9BXJ1](#)
NCBI
[Q9BXJ1](#)
Dilution Range

This recombinant protein can be used for biological assays. For research use only.

Expiration Date

6 months from date of receipt.