


Product Datasheet

PRMT7 antibody (orb12297)

Description	Goat polyclonal to PRMT7.
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRMT7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms NP_061896.1;NP_001171753.1.
Protein Sequence	PRFGEINDQDRTDR
Clonality	Polyclonal
MW	78.5; 73.2
Entrez	54496
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using PRMT7 antibody (orb12297).