

Product Datasheet

CXCL14 Recombinant Protein (orb1229817)

Description

CXCL14 Recombinant Protein

Target

CXCL14

Preservatives

Lyophilized from a 0.2 um filtered solution of 20mM TrisHCl, 1M NaCl, pH 8.5. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH2O.

Form/Appearance

Lyophilized

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at 4-7°C for 2-7 days.
Aliquots of reconstituted samples are stable at -20°C for 3 months.

Tag

Tag Free

Note

For research use only

Application notes

This recombinant protein can be used for biological assays. For research use only.

Protein Sequence

Ser35-Glu111

Purity

Greater than 95% as determined by reducing SDS-PAGE.
Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.

Source

E. coli

MW

9.4 kD

Uniprot ID
O95715
NCBI
O95715
Dilution Range

This recombinant protein can be used for biological assays. For research use only.

Expiration Date

6 months from date of receipt.