

Product Datasheet

Tim-4 Recombinant Protein (orb1229866)

Description	Tim-4 Recombinant Protein
Target	Timd4
Preservatives	Lyophilized from 0.2um-filtered solution in PBS. Reconstitute with 100 µl sterile water. Add 1X PBS to the desired protein concentration.
Form/Appearance	Lyophilized
Storage	Stable for at least 1 year after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Tag	Fc Tag
Note	For research use only
Application notes	This recombinant proteins is for research use only.
Protein Sequence	The extracellular domain of mouse Tim-4 (aa 22-279) is fused to the N-terminus of the Fc region of human IgG1.
Purity	≥98% (SDS-PAGE). Endotoxin level is less than 0.06EU/ µg protein (LAL test; Lonza).
Source	CHO cells
Uniprot ID	O53649
NCBI	NP_848874
Dilution Range	This recombinant proteins is for research use only.
Expiration Date	6 months from date of receipt.