

Product Datasheet

Hepatitis E virus ORF 3 Recombinant Protein (orb1230308)

Description	Hepatitis E virus ORF 3 Recombinant Protein
Tested Applications	ELISA
Preservatives	20 mM Tris, pH 8.0 + 8M urea + 10 mM BME
Form/Appearance	liquid
Concentration	batch dependent
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Tag	beta-galactosidase
Note	For research use only
Application notes	This recombinant protein can be used ELISA. For research use only.
Protein Sequence	amino acids 92 - 123
Purity	~ 95% by SDS-PAGE
Source	E. coli
Dilution Range	This recombinant protein can be used ELISA. For research use only.
Expiration Date	6 months from date of receipt.