

Product Datasheet

Hepatitis C virus NS5a Recombinant Protein (orb1230310)

Description	Hepatitis C virus NS5a Recombinant Protein
Tested Applications	ELISA
Preservatives	50 mM Tris, pH 8.0 + EDTA
Form/Appearance	liquid
Concentration	batch dependent
Storage	Store in working aliquots at 2-8°C. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Tag	GST
Note	For research use only
Application notes	This recombinant protein can be used ELISA. For research use only.
Protein Sequence	NS5a region
Purity	~ 90% by SDS-PAGE
Source	E. coli
MW	38 kDa
Dilution Range	This recombinant protein can be used ELISA. For research use only.
Expiration Date	6 months from date of receipt.