

Product Datasheet

Hepatitis C virus E2 Recombinant Protein (orb1230313)

Description	Hepatitis C virus E2 Recombinant Protein
Tested Applications	ELISA
Preservatives	25 mM Tris-HCl, 0.1% SDS, pH 8.3
Form/Appearance	liquid
Concentration	batch dependent
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Tag	none
Note	For research use only
Application notes	This recombinant protein can be used ELISA. For research use only.
Protein Sequence	NA
Purity	~ 95% by SDS-PAGE
Source	E. coli
MW	28 kDa
Dilution Range	This recombinant protein can be used ELISA. For research use only.
Expiration Date	6 months from date of receipt.