

Product Datasheet

Hepatitis C virus NS3 Recombinant Protein (orb1230316)

Description	Hepatitis C virus NS3 Recombinant Protein
Tested Applications	ELISA
Preservatives	25 mM Tris, pH 8.3 + 0.1% SDS
Form/Appearance	liquid
Concentration	batch dependent
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Tag	none
Note	For research use only
Application notes	This recombinant protein can be used ELISA. For research use only.
Protein Sequence	amino acids 1230 - 1501
Purity	~ 90% by SDS-PAGE
Source	E. coli
MW	42 kDa
Dilution Range	This recombinant protein can be used ELISA. For research use only.
Expiration Date	6 months from date of receipt.