

Product Datasheet

Hepatitis C virus NS4 Recombinant Protein (orb1230317)

Description	Hepatitis C virus NS4 Recombinant Protein
Tested Applications	ELISA, WB
Preservatives	10 mM Na ₂ CO ₃ /10 mM EDTA/14 mM beta-ME/0.05% tween 20
Form/Appearance	liquid
Concentration	batch dependent
Storage	Store in working aliquots at 2-8°C. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Tag	beta-galactosidase (114 kDa)
Note	For research use only
Application notes	This recombinant protein can be used for WB and ELISA. For research use only.
Protein Sequence	NS4a+b; amino acids 1658 -1863 (19 kDa)
Purity	~ 95% by SDS-PAGE
Source	E. coli
MW	133 kDa
Dilution Range	This recombinant protein can be used for WB and ELISA. For research use only.
Expiration Date	6 months from date of receipt.