

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Recombinant protein (orb1231007)

Description	SARS-CoV-2 (COVID-19) Nucleocapsid Recombinant protein
Preservatives	CoV 2019 protein solution is supplied in PBS and 10mM K ₂ CO ₃ .
Form/Appearance	Liquid
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles.
Tag	His Tag
Note	For research use only
Protein Sequence	The E.Coli derived recombinant protein contains the Coronavirus 2019 full-length nucleoprotein S- gene: Gene bank- MN908947, fused to 6xHis tag at C-terminal and having a Mw. Of 48 kDa as appears on SDS-PAGE.
Purity	> 95% pure as determined SDS- PAGE.
Source	E. coli
Expiration Date	6 months from date of receipt.