

Product Datasheet

SARS-CoV-2 (COVID-19) Spike Glycoprotein-S1, Recombinant protein (orb1231009)

Description	SARS-CoV-2 (COVID-19) Spike Glycoprotein-S1, Recombinant protein
Preservatives	nCoV-S1 protein solution is supplied in DPBS.
Form/Appearance	Liquid
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles.
Tag	Fc tag
Note	For research use only
Protein Sequence	The HEK293 derived recombinant protein contains the Novel Coronavirus 2019-nCoV Spike Glycoprotein S1, Wuhan-Hu-1 strain, amino acids 1-674 fused to Fc tag at C-terminal.
Purity	> 85% pure as determined SDS-PAGE. Purified by Protein-G chromatographic technique.
Source	HEK293 cells
Expiration Date	6 months from date of receipt.