


Product Datasheet

SARS-CoV-2 (COVID-19) Delta Variant (B.1.617.2) Spike RBD (L452R, T478K) Recombinant Protein (orb1231024)

Description	SARS-CoV-2 (COVID-19) Delta Variant (B.1.617.2) Spike RBD (L452R, T478K) Recombinant Protein
Preservatives	Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH ₂ O to the concentration with 1mg/ml.
Form/Appearance	Lyophilized
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
Tag	His Tag
Note	For research use only
Application notes	Immunogen
Protein Sequence	RBD with mutations as L452R and T478K.
Purity	> 90% as determined by SDS-PAGE
Source	Mammalian Cells
NCBI	YP_009724390.1, L452, T478
Dilution Range	Immunogen
Expiration Date	6 months from date of receipt.



SDS-PAGE Image