



Product Datasheet


**SARS-CoV-2 (COVID-19) Full-Length Eta Variant
(B.1.525, UK) Spike Recombinant Protein
(orb1231034)**

Description

SARS-CoV-2 (COVID-19) Full-Length Eta



A: SDS-
Page of
purified
SPIKE in
deterge...



	Protein
Tested Applications	ELISA
Target	S
Preservatives	20 mM Hepes pH 7.5; 150 mM NaCl, 0.001 % LMNG
Form/Appearance	Liquid
Concentration	batch dependent
Storage	-80°C. Avoid freeze-thaw cycles
Tag	C-terminal Rho1D4
Note	For research use only
Application notes	This recombinant protein can be used for ELISA assays, Ligand binding assays (e.g. SPR), Biochemical and biophysical analyses. For research use only.
Protein Sequence	Full-length sequence (aa 1 - 1273), del 69-70; del 144; Q52R; E484K; Q677H; F888Lfurin cleavage site "RRAR" mutated to "GSAG"; K986P; V987PC-terminal Rho1D4 tag fused with spacer "GSSG" to protein sequence
Purity	> 98% as determined by SDS-PAGE, see Fig. 1 A and B. PureCube Rho1D4 Agarose.
Source	HEK293 cells
MW	1286 amino acids (including Rho1D4 tag and linker); 142114 Da
Uniprot ID	P0DTC2
Dilution Range	This recombinant protein can be used for ELISA assays, Ligand binding assays (e.g. SPR), Biochemical and biophysical analyses. For research use only.
Expiration Date	6 months from date of receipt.