


Product Datasheet

**SARS-CoV-2 (COVID-19) Gamma Variant (P.1,
Brazil) Spike S1 (RBD) Variant Recombinant
Protein (orb1231045)**

Description

SARS-CoV-2 (COVID-19) Gamma Variant



SARS-CoV-2 Spike Protein S1 (RBD):Fc (hu...)

	Recombinant Protein
Target	S
Preservatives	Lyophilized. Contains PBS.
Form/Appearance	Lyophilized
Concentration	1mg/ml after reconstitution.
Storage	Short Term Storage: +4°C, Long Term Storage: -20°C, Handling Advice: After opening, prepare aliquots and store at -20°C. Avoid freeze/thaw cycles. Centrifuge lyophilized vial before opening and reconstitution. For maximum product recovery after thawing, centrifuge the vial before opening the cap. Use/Stability: Stable for at least 6 months after receipt when stored at -20°C. Working aliquots are stable for up to 3 months when stored at -20°C.
Note	For research use only
Protein Sequence	Receptor-binding domain (RBD) of SARS-CoV-2 Spike protein S1 (aa 319-541) containing the mutations K417T, E484K, N501Y is fused to the N-terminus of the Fc region of human IgG1.
Purity	≥95% (SDS-PAGE)
Source	HEK293 cells
MW	~60kDa (SDS-PAGE); ~115 kDa (Dimer) (SEC)
Uniprot ID	P0DTC2
NCBI	QHD43416.1
Expiration Date	6 months from date of receipt.