


Product Datasheet

CHRM2 antibody (orb12313)

Description	Goat polyclonal to CHRM2.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CHRM2 (aa352-364)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Mouse Brain lysates (calculated MW of 51.5kDa according to Mouse NP_987076.2). Recommended concentration: 0.5-2µg/ml.
Protein Sequence	QNGDEKQNIVARK
Clonality	Polyclonal
MW	51.7
Entrez	1129
Dilution Range	ELISA: 1:128000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Brain lys...