


Product Datasheet

SARS-CoV-2 (COVID-19) S1 Recombinant Protein NTD (orb1231454)

Description

SARS-CoV-2 (COVID-19) S1 Recombinant Protein



Tested Applications

ELISA, WB

Target

S

Preservatives

PBS, pH 7.3

Form/Appearance

Lyophilized

Storage

Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.

Tag

His Tag

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 34.9 kDa. The protein migrates as 45-65 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.

Protein Sequence

Ser 13 - Leu 303

Purity

> 90% as determined by SDS-PAGE.

Source

HEK293 cells

MW

34.9 kDa

Uniprot ID

[P0DTC2](#)

NCBI

[QHD43416.1](#)


Dilution Range

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 34.9 kDa. The protein migrates as 45-65 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.

Expiration Date

6 months from date of receipt.

SARS-CoV-2
S1 protein
NTD, His
Tag on
SD...



Immobilized
SARS-CoV-2
S1 protein
NTD, H...