



Product Datasheet

RIF1 antibody (orb12315)


Description	Goat polyclonal to RIF1.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	RIF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms NP_060621.3;NP_001171134.1. The immunizing peptide sequence has been optimised for cross reactivity in mouse. Reported variants represent identical protein: NP_001171135.1, NP_001171136.1, NP_001171134.1
Protein Sequence	RNLQSRWRSPSHEN
Clonality	Polyclonal
MW	274; 272
Entrez	55183
Dilution Range	ELISA: 1:128000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa and Mouse ...



Immunofluorescence analysis of NIH3T3 ce...



Immunofluorescence analysis of U2OS cell...