

## Product Datasheet

### Human ACE2 Recombinant Protein (orb1231569)

|                          |                                                                                                                                                                         |                                 |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| <b>Description</b>       | Human ACE2 Recombinant Protein                                                                                                                                          | Image not found or type unknown |
| <b>Preservatives</b>     | Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH2O to the concentration with 1mg/ml.                                                  | SDS-PAGE<br>Image               |
| <b>Form/Appearance</b>   | Lyophilized                                                                                                                                                             |                                 |
| <b>Storage</b>           | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt. |                                 |
| <b>Tag</b>               | C-His Tag                                                                                                                                                               |                                 |
| <b>Note</b>              | For research use only                                                                                                                                                   |                                 |
| <b>Application notes</b> | Immunogen                                                                                                                                                               |                                 |
| <b>Protein Sequence</b>  | Ser19-Asp615                                                                                                                                                            |                                 |
| <b>Purity</b>            | > 90% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining.                                                                                   |                                 |
| <b>Source</b>            | Mammalian Cells                                                                                                                                                         |                                 |
| <b>MW</b>                | 71.54kDa                                                                                                                                                                |                                 |
| <b>NCBI</b>              | <a href="#">NP_001358344.1</a> , <a href="#">Q9BYF1</a>                                                                                                                 |                                 |
| <b>Dilution Range</b>    | Immunogen                                                                                                                                                               |                                 |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                          |                                 |