

Product Datasheet

SARS-CoV-2 (COVID-19) S1 Recombinant Protein (orb1231580)

Description	SARS-CoV-2 (COVID-19) S1 Recombinant Protein
Preservatives	Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH2O to the concentration with 1mg/ml.
Form/Appearance	Lyophilized
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
Tag	C-His Tag
Note	For research use only
Application notes	Immunogen
Protein Sequence	Met1-Tyr695
Purity	> 80% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining.
Source	HEK293 Cells
MW	77.23kDa
NCBI	QHD43416.1
Dilution Range	Immunogen
Expiration Date	6 months from date of receipt.