

Product Datasheet

SARS-CoV-2 (COVID-19) Spike RBD Recombinant Protein (orb1231588)

Description

SARS-CoV-2 (COVID-19) Spike RBD

Tested Applications

ELISA, WB

Preservatives

PBS

Form/Appearance

Liquid

Storage

The protein is stable in a liquid state at -70°C for 12 months. Avoid repeated freeze-thaw cycles.

Tag

His tag at C-terminus

Note

For research use only

Protein Sequence

Arg319-Phe541

Purity

95%

Source


HEK293 cells

MW


24.7 kD

Expiration Date


6 months from date of receipt.



Binding Assay of ACE2 and SARS-CoV-2 (COV...



SDS-PAGE Validation with SARS-CoV-2 (COV...



Western Blot Validation with SARS-CoV-2 ...