


Product Datasheet

SARS-CoV-2 (COVID-19) Papain-like Protease (orb1231605)

Description	SARS-CoV-2 (COVID-19) Papain-like Protease
Tested Applications	WB
Preservatives	20mM Tris, 20% Glycerol, 10mM β-Me, pH 7.5
Form/Appearance	Liquid
Storage	This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.
Tag	His tag at the N-terminus
Note	For research use only
Application notes	This recombinant protein can be used for SDS-PAGE.
Protein Sequence	Glu1564-Lys1878
Purity	> 90% by SDS-PAGE.
Source	E. coli
MW	36.6 kDa
NCBI	YP_009725299.1
Dilution Range	This recombinant protein can be used for SDS-PAGE.
Expiration Date	6 months from date of receipt.



SARS-CoV-2
(COVID-19)
papain-like
protea...