



Product Datasheet

SARS-CoV-2 (COVID-19) Spike S1 Recombinant Protein (orb1231606)

Description	SARS-CoV-2 (COVID-19) Spike S1 Recombinant Protein
Tested Applications	ELISA, WB
Target	S
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week.
Tag	His tag at the C-terminus
Note	For research use only
Application notes	This recombinant protein can be used for SDS-PAGE.
Protein Sequence	Val11-Arg683
Purity	> 90% by SDS-PAGE.
Source	HEK293 cells
MW	76 kDa
Uniprot ID	P0DTC2
NCBI	YP_009724390.1
Dilution Range	This recombinant protein can be used for SDS-PAGE.
Expiration Date	6 months from date of receipt.



Immobilized SARS-CoV-2 (COVID-19) spike ...



SARS-CoV-2 (COVID-19) spike S1 recombina...