


Product Datasheet

SARS-CoV-2 (COVID-19) 3C-like Proteinase (orb1231608)

Description	SARS-CoV-2 (COVID-19) 3C-like Proteinase
Tested Applications	WB
Target	ORF1ab
Preservatives	20mM HEPES, 120mM NaCl, 10% Glycerol, 2mM TCEP, pH 7.4
Form/Appearance	Liquid
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles.
Tag	His, Avi tag at the N-terminus
Note	For research use only
Application notes	This recombinant protein can be used for SDS-PAGE.
Protein Sequence	Ser1-Gln306
Purity	> 95% by SDS-PAGE.
Source	E. coli
MW	36.9 kDa
Uniprot ID	P0DTD1
NCBI	YP_009725301.1
Dilution Range	This recombinant protein can be used for SDS-PAGE.
Expiration Date	6 months from date of receipt.



SARS-CoV-2
(COVID-19)
3C-like
proteinase...