



Product Datasheet

Human ACE2 Recombinant Protein (orb1231610)

Description	Human ACE2 Recombinant Protein
Tested Applications	ELISA, WB
Target	ACE2
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.
Tag	His, Avi tag at the C-terminus
Note	For research use only
Application notes	This recombinant protein can be used for SDS-PAGE.
Protein Sequence	Gln18- Ser740
Purity	> 95% by SDS-PAGE.
Source	HEK293 cells
MW	86.2 kDa
Uniprot ID	Q9BYF1
NCBI	Q9BYF1
Dilution Range	This recombinant protein can be used for SDS-PAGE.
Expiration Date	6 months from date of receipt.



Immobilized
human
ACE2
recombinant
prote...



Human
ACE2
recombinant
protein was
deter...