

Product Datasheet

PEX10 (CT) Peptide (orb1231632)

Description	PEX10 (CT) Peptide
Tested Applications	Blocking
Target	PEX10
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store PEX10°CT) Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	PEX10 (CT) Peptideis used for blocking the activity of thePEX10 (CT) Antibody
Uniprot ID	O60683
NCBI	O60683
Expiration Date	6 months from date of receipt.